

# BIDS Curation with fw-heudiconv

Comprehensive guide to curating your imaging data on Flywheel

Tinashe Michael Tapera

August 2019


# To do:



1. Introduction to `fw-heudiconv`
  2. How does it work?
  3. Hands-on walkthrough
  4. Workshop (optional)
-

# Introduction to fw-heudiconv

1. Flywheel & BIDS
2. heudiconv + Flywheel = **fw-heudiconv**



# Flywheel & BIDS

- **Flywheel:** *Next generation informatics platform for biomedical research and collaboration*
- Eliminates challenges of traditional neuroscience research:
  - Boring
  - Difficult
  - Error-prone
  - Not directly related to neuroscience (e.g. submitting compute jobs, compiling software)
- **BIDS:** *A simple and intuitive way to organize and describe your neuroimaging and behavioral data*
- Enhances traditional neuroscience research:
  - Easy collaboration between labs and databases
  - Easy consumption by conversion & analysis software applications

# Flywheel & BIDS: “Traditional” Research Model

1. Identify researcher
2. Download scans from scanner
3. Organize your data
4. Choose analysis software
5. Write custom processing pipeline code
6. Conduct quality control
7. Pass off analyses to next stakeholder

Trainees have different skill levels

Conversion software varies, and choices matter

Directory structures vary across labs

Multiple dependencies of software versions

Code poorly documented, no version control.

Quality control often idiosyncratic

Documentation usually poor, does not keep pace with versions.

# Flywheel & BIDS: Reproducible Neuroscience Model

1. Data pulled straight from scanner (*Flywheel reaper*)
2. Data curation with BIDS that describes data and acquisition parameters (**fw-heudiconv**)
3. High-performing, benchmarked pre-processing pipelines that auto-configure to data acquisition (*BIDS-apps\**)
4. Automated, reproducible quality control
5. Guaranteed reproducibility and “glass box” code with wrapped dependencies (*Docker containers*)
6. Ease of sharing across labs and projects (*Flywheel gears*)
7. Unlimited scalability (*Google cloud platform*)

\* Demos in the next meeting



So what does data curated in BIDS look like...?

# Flywheel & BIDS: BIDS-ified Data

Key features:

- `_{key}-{value}_` pairs separated by underscores
- Entity hierarchy: Project > Subject > Session > Acquisition > File
- JSON sidecars for parseable metadata

Using BIDS, the data describes itself!

**sub**-{value}, **ses**-{value}, **acq**-{value}, **ce**-{value},  
**rec**-{value}, **dir**-{value}, **run**-{index}, **mod**-{value}...

**fw-heudiconv** helps us accomplish this on  
Flywheel!

```
.bidsignore
dataset_description.json
sub-106246/
  ses-1/
    anat/
      sub-106246_ses-1_T1w.json
      sub-106246_ses-1_T1w.nii.gz
    fmap/
      sub-106246_ses-1_magnitude1.json
      sub-106246_ses-1_magnitude1.nii.gz
      sub-106246_ses-1_magnitude2.json
      sub-106246_ses-1_magnitude2.nii.gz
    func/
      sub-106246_ses-1_task-frac2back.json
      sub-106246_ses-1_task-frac2back.nii.gz
      sub-106246_ses-1_task-idemo.json
      sub-106246_ses-1_task-idemo.nii.gz
```



# How does it work?


---

1. Flywheel SDK and Gears
2. Writing heuristics:
  - a. `create_key()`
  - b. `infotodict()`
  - c. `MetadataExtras`
  - d. `IntendedFor`
  - e. `ReplaceSubject()` &  
`ReplaceSession()`

# SDK & Gears

- **Software Development Kit:** *a set of tools, libraries, relevant documentation, code samples, processes, and guides that allow developers to create software applications on some platform*
  - Available in Python, Matlab, R
- **Gears:** *ready-to-use pre-processing or analysis pipeline applications; “an algorithm that has been packaged in a way that allows it to be run and managed within the Flywheel Platform”*
  - Containerised (through Docker), so almost anything is possible!
  - Metadata is preserved
  - Code is visible
  - Version controlled

# SDK & Gears





So when are we going to see some code already?



HeuDiConv

+

***FLYVHEEL***

“a flexible DICOM converter for  
organizing brain imaging data into  
structured BIDS directory layouts”

=

**BBL**  
Brain Behavior Lab

fw-heudiconv

# Heuristics with `fw-heudiconv`

- `fw-heudiconv` is written in Python, and builds upon `heudiconv` principles (i.e. heuristics)
- **Heuristic**: *a discrete set of rules that differentiates imaging files by their DICOM header information*
- Like the original `heudiconv`, `fw-heudiconv` looks through your scans on Flywheel and renames them based on the ruleset you specify using **boolean logic** in your heuristic file.

Let's look at some Python functions that make up a heuristic file...

# Basic Functions: `create_key()`

This function defines template keys for the different scan types and the naming convention they will take. `fw-heudiconv` uses string formatting to insert the values.

```
tlw = create_key('sub-{subject}/{session}/anat/sub-{subject}_{session}_T1w')
```

You can create as many keys as you want for your scans (the more strict/granular the better)

```
nback_HiConHiLoWMgated_run1 = create_key(
```

```
'sub-{subject}/{session}/func/sub-{subject}_{session}_'
```

```
'task-nback_acq-HiConHiLoWMgated_run-01_bold')
```

# Basic Functions: `infotodict()`

This function loops over the scans and their metadata, and assigns each scan a template key in a variable. Using a python dictionary like this:

```
info = {t1w:[], nback_HiConHiLoWMgated_run1:[]}
```

`fw-heudiconv` will add a scan to the correct list if it meets some logical criteria:

```
if "anat_t1w" in protocol:
```

```
    info[t1w].append(s.series_id)
```

```
elif "HiConHiLoWMgated" in s.protocol_name and "M" in s.image_type:
```

```
    info[nback_HiConHiLoWMgated_run1].append(s.series_id)
```



# Extensions: **MetadataExtras()**

**fw-heudiconv** functionality is easily extensible for Flywheel specific scenarios.

**MetadataExtras()** is a variable that allows users to hard code metadata to key templates:

```
MetadataExtras = {  
  
    nback_HiConHiLoWMgated_run1: { "TaskName": "n-back" }  
  
}
```

These metadata show up in the JSON sidecar file on export, and can overwrite pre-existing Flywheel data or create new fields

# Extensions: **IntendedFor()**

Using fieldmaps? BIDS-enabled applications such as ***fMRIPrep*** need to be told what scans these fieldmaps are correcting for. This is easy to do with is a variable that specifies these in **fw-heudiconv**.

```
IntendedFor = {
```

```
    fmap: [
```

```
        '{session}/func/sub-{subject}_{session}_task-nback_acq-HiConHiLoWMgated_run-01_bold'
```

```
    ]
```

```
}
```

**fw-heudiconv** will check to make sure that each of these files exists post-hoc

## Extensions: **Replace\*** ()

What if you have subject or session labels that you want to correct in BIDS?

**fw-heudiconv** has a solution for this too. **fw-heudiconv** has dedicated functionality for parsing and replacing subject labels or session labels, so you can write a custom function in the heuristic that tells **fw-heudiconv** how to handle each case:

```
def ReplaceSubject(subj_label):
```

```
    return str(int(subj_label))
```

```
def ReplaceSession(sess_label):
```

```
    return str(int(sess_label + 1))
```

# Walkthrough

---

1. Querying a project
2. Curating your data
3. Exporting your data

# Walkthrough: Querying


First step to curating your data: *understanding your data!*

Use the tabulation tool to get a tabular view of the unique dicom sequences in your dataset:

```
$ fw-heudiconv-tabulate --project gear_testing --path <where to download the table>  
--subject <optional list of subjects> --session <optional list of sessions> --dry_run <to  
print and not download>
```

Unique fields include protocol names, sequence names, TR, TE, motion correction, and derived images.

# Walkthrough: Querying



```
fw-heudiconv-demo — ttapera@chead:~ — -bash — 74x22
(flywheel) Tinashes-Air:fw-heudiconv-demo mac$ time fw-heudiconv-tabulate
--project gear_testing --path . --subject exampleSubject
/Users/mac/anaconda2/envs/flywheel/lib/python3.7/site-packages/fw_heudicon
v/query.py:4: UserWarning: The DICOM readers are highly experimental, unst
able, and only work for Siemens time-series at the moment
Please use with caution. We would be grateful for your help in improving
them
    from nibabel.nicom.dicomwrappers import wrapper_from_data
INFO: Querying Flywheel server...
INFO: Done!

real    0m21.267s
user    0m3.279s
sys     0m1.103s
(flywheel) Tinashes-Air:fw-heudiconv-demo mac$
```

# Walkthrough: Querying

[illegible]

# Walkthrough: Curating

Next we begin curating the data into BIDS.

1. Add the basic heuristic functions to a new python file.
2. Use the curation tool in *dry\_run* mode to test out and see what changes would be applied if you used the heuristic:

```
$ fw-heudiconv-curate --project gear_testing --heuristic <path to heuristic file> --subject  
<optional list of subjects> --session <optional list of sessions> --dry_run <to print and  
not apply changes>
```



# Walkthrough: Curating

```
fw-heudiconv-demo -- ttapera@chead:~ -- bash -- 80x56

(flywheel) Tinashes-Air:fw-heudiconv-demo mac$ fw-heudiconv-curate --project gear_testing --heuristic demo_heuristic.py --dry_run --subject exampleSubject /Users/mac/anaconda2/envs/flywheel/lib/python3.7/site-packages/fw-heudiconv/queries.py:4: UserWarning: The DICOM wrappers are highly experimental, unstable, and on ly work for Siemens time-series at the moment
Please use with caution. We would be grateful for your help in improving them
from nibabel.nicom.dicomwrappers import wrapper_from_data
INFO: Querying Flywheel server...
DEBUG: Found project: gear_testing (5c8937fdd93e3802e025e2b)
DEBUG: Found sessions:
      exampleSession (5d488062c3758a0043d51157)
      exampleSession2 (5d488060c3758a0043d51155)
DEBUG: Found SeqInfos:
MPRAGE_T1100_ipat2:
      [TR=1.81 TE=0.00351 shape=(256, 192, 160, -1) image_type=('ORIGINAL', 'PRIMARY', 'M', 'ND', 'NORM')] (5d488064c3758a0047d684bc)

ep2d_se_pcasl_PHC_1200ms_MoCo:
      [TR=4 TE=0.829 shape=(480, 480, 80, -1) image_type=('ORIGINAL', 'PRIMARY', 'M', 'ND', 'MOCO', 'MOSAIC')] (5d488063c3758a0049d7657f)

bbl1_cardB0_178:
      [TR=2 TE=0.825 shape=(384, 384, 178, -1) image_type=('ORIGINAL', 'PRIMARY', 'M', 'ND', 'MOSAIC')] (5d488063c3758a003ed47edb)

bbl1_faceA0_178:
      [TR=2 TE=0.825 shape=(384, 384, 178, -1) image_type=('ORIGINAL', 'PRIMARY', 'M', 'ND', 'MOSAIC')] (5d488063c3758a0047d684b9)

bbl1_faceB0_178:
      [TR=2 TE=0.825 shape=(384, 384, 178, -1) image_type=('ORIGINAL', 'PRIMARY', 'M', 'ND', 'MOSAIC')] (5d488063c3758a0045d59f32)

ep2d_bold_MGH:
      [TR=3 TE=0.83 shape=(504, 504, 124, -1) image_type=('ORIGINAL', 'PRIMARY', 'M', 'ND', 'MOSAIC')] (5d488063c3758a0045d59f30)

ep2d_se_pcasl_PHC_1200ms:
      [TR=4 TE=0.829 shape=(480, 480, 80, -1) image_type=('ORIGINAL', 'PRIMARY', 'M', 'ND', 'MOSAIC')] (5d488063c3758a0049d76582)

epi_singlerep:
      [TR=2 TE=0.83 shape=(384, 384, 1, -1) image_type=('ORIGINAL', 'PRIMARY', 'M', 'ND', 'MOSAIC')] (5d488063c3758a003dd47a94)

localizer_32channel:
      [TR=0.0086 TE=0.004 shape=(512, 512, 3, -1) image_type=('ORIGINAL', 'PRIMARY', 'M', 'ND', 'NORM')] (5d488063c3758a0049d76581)

localizer_32channel:
      [TR=0.0086 TE=0.004 shape=(512, 512, 3, -1) image_type=('ORIGINAL', 'PRIMARY', 'M', 'ND', 'NORM')] (5d488063c3758a0049d76581)

localizer_32channel:
      [TR=0.0086 TE=0.004 shape=(512, 512, 3, -1) image_type=('ORIGINAL
```

```
fw-heudiconv-demo -- ttapera@chead:~ -- bash -- 80x56

', 'PRIMARY', 'M', 'MB', 'ND', 'MOSAIC')] (5d488062c3758a0047d684b6)

ep2d_single:
      [TR=3 TE=0.83 shape=(448, 448, 1, -1) image_type=('ORIGINAL', 'PRIMARY', 'M', 'ND', 'MOSAIC')] (5d488062c3758a0049d76574)

localizer:
      [TR=0.02 TE=0.005 shape=(192, 192, 3, -1) image_type=('ORIGINAL', 'PRIMARY', 'M', 'ND')] (5d488063c3758a0049d7657b)

localizer:
      [TR=0.02 TE=0.005 shape=(192, 192, 3, -1) image_type=('ORIGINAL', 'PRIMARY', 'M', 'ND')] (5d488063c3758a0049d7657b)

localizer:
      [TR=0.02 TE=0.005 shape=(192, 192, 3, -1) image_type=('ORIGINAL', 'PRIMARY', 'M', 'ND')] (5d488063c3758a0049d7657b)

INFO: Loading heuristic file...
INFO: Applying heuristic to query results...
Series not recognized!: ep2d_se_pcasl_PHC_1200ms_MoCo ep2d_se_pcasl_PHC_1200ms_MoCo_4.nii.gz
Series not recognized!: bbl1_cardB0_178 bbl1_cardB0_178_9.nii.gz
Series not recognized!: bbl1_faceA0_178 bbl1_faceA0_178_8.nii.gz
Series not recognized!: bbl1_faceB0_178 bbl1_faceB0_178_10.nii.gz
Series not recognized!: ep2d_bold_MGH ep2d_bold_MGH_13.nii.gz
Series not recognized!: ep2d_se_pcasl_PHC_1200ms ep2d_se_pcasl_PHC_1200ms_3.nii.gz
Series not recognized!: epi_singlerep epi_singlerep_6.nii.gz
Series not recognized!: localizer_32channel localizer_32channel_1_100003.nii.gz
Series not recognized!: localizer_32channel localizer_32channel_1_100002.nii.gz
Series not recognized!: localizer_32channel localizer_32channel_1_100001.nii.gz
Series not recognized!: ep2d_effort1_236 ep2d_effort1_236_9.nii.gz
Series not recognized!: ep2d_effort2_236 ep2d_effort2_236_18.nii.gz
Series not recognized!: ep2d_effort3_1416 ep2d_effort3_1416_11.nii.gz
Series not recognized!: ep2d_single ep2d_single_6.nii.gz
Series not recognized!: localizer_localizer_1_100003.nii.gz
Series not recognized!: localizer_localizer_1_100002.nii.gz
Series not recognized!: localizer_localizer_1_100001.nii.gz
DEBUG:
MPRAGE_T1110 ipat2_moco3_4.nii.gz
      sub-exampleSubject ses-exampleSession2_T1w.nii.gz -> sub-exampleSubject/ses-exampleSession2/anat/sub-exampleSubject_ses-exampleSession2_T1w.nii.gz
DEBUG:
MPRAGE_T1100 ipat2_2.nii.gz
      sub-exampleSubject ses-exampleSession_T1w.nii.gz -> sub-exampleSubject/ses-exampleSession/anat/sub-exampleSubject_ses-exampleSession_T1w.nii.gz
DEBUG:
MPRAGE_NAVprotocol_2.nii.gz
      sub-exampleSubject ses-exampleSession2_T1w.nii.gz -> sub-exampleSubject/ses-exampleSession2/anat/sub-exampleSubject_ses-exampleSession2_T1w.nii.gz
DEBUG:
MPRAGE_T1110 ipat2_moco3_3.nii.gz
      sub-exampleSubject ses-exampleSession2_T1w.nii.gz -> sub-exampleSubject/ses-exampleSession2/anat/sub-exampleSubject_ses-exampleSession2_T1w.nii.gz
(flywheel) Tinashes-Air:fw-heudiconv-demo mac$ []
```

# Walkthrough: Curating

3. Add more keys to capture the sequences not recognised by your logic
4. Rinse and repeat
5. Add hardcoded metadata, fieldmap intentions, and subject/session label replacements as needed

# Walkthrough: Curating

```
fw-heudiconv-demo — ttapera@chead:~ — -bash — 172x66

Series not recognized!: ep2d_effort1_236 ep2d_effort1_236_9.nii.gz
Series not recognized!: ep2d_effort2_236 ep2d_effort2_236_10.nii.gz
Series not recognized!: ep2d_effort3_1416 ep2d_effort3_1416_11.nii.gz
Series not recognized!: ep2d_single ep2d_single_6.nii.gz
Series not recognized!: localizer_localizer_1_i00003.nii.gz
Series not recognized!: localizer_localizer_1_i00002.nii.gz
Series not recognized!: localizer_localizer_1_i00001.nii.gz
INFO: Processing IntendedFor fields based on heuristic file
DEBUG: Intention map: defaultdict(<class 'list'>, {(('sub-(subject)/(session)/fmap/sub-(subject)/(session)_magnitude(item)', ('nii.gz',)), None): [('session)/func/sub-(subject)_task-rest_bold.nii.gz', ('session)/func/sub-(subject)_task-face_run-01_bold.nii.gz', ('session)/func/sub-(subject)_task-face_run-02_bold.nii.gz', ('session)/func/sub-(subject)_task-card_run-01_bold.nii.gz', ('session)/func/sub-(subject)_task-card_run-02_bold.nii.gz'], ('sub-(subject)/(session)/fmap/sub-(subject)_phase(item)', ('nii.gz',)), None): [('session)/func/sub-(subject)_task-rest_bold.nii.gz', ('session)/func/sub-(subject)_task-face_run-01_bold.nii.gz', ('session)/func/sub-(subject)_task-face_run-02_bold.nii.gz', ('session)/func/sub-(subject)_task-card_run-01_bold.nii.gz', ('session)/func/sub-(subject)_task-card_run-02_bold.nii.gz']})
INFO: Processing Metadata fields based on heuristic file
DEBUG: Metadata extras: defaultdict(<class 'list'>, {(('sub-(subject)/(session)/func/sub-(subject)_task-rest_bold', ('nii.gz',)), None): {'TaskName': 'rest'}, ('sub-(subject)/(session)/func/sub-(subject)_task-card_run-01_bold', ('nii.gz',)), None): {'TaskName': 'card'}, ('sub-(subject)/(session)/func/sub-(subject)_task-card_run-02_bold', ('nii.gz',)), None): {'TaskName': 'card'}, ('sub-(subject)/(session)/func/sub-(subject)_task-face_run-01_bold', ('nii.gz',)), None): {'TaskName': 'face'}, ('sub-(subject)/(session)/func/sub-(subject)_task-face_run-02_bold', ('nii.gz',)), None): {'TaskName': 'face'}, ('sub-(subject)/(session)/asl/sub-(subject)_asl', ('nii.gz',)), None): {'ASLContext': 'Label_Control', 'LabelingType': 'PCASL'}})
DEBUG:
MPRAGE_T11100_ipat2_2.nii.gz
sub-01_ses-1_T1w.nii.gz -> sub-01/ses-1/anat/sub-01_ses-1_T1w.nii.gz
DEBUG:
MPRAGE_T11110_ipat2_moco3_3.nii.gz
sub-01_ses-2_T1w.nii.gz -> sub-01/ses-2/anat/sub-01_ses-2_T1w.nii.gz
DEBUG:
MPRAGE_NAVprotocol_2.nii.gz
sub-01_ses-2_T1w.nii.gz -> sub-01/ses-2/anat/sub-01_ses-2_T1w.nii.gz
DEBUG:
MPRAGE_T11110_ipat2_moco3_4.nii.gz
sub-01_ses-2_T1w.nii.gz -> sub-01/ses-2/anat/sub-01_ses-2_T1w.nii.gz
DEBUG:
ep2d_se_pcasl_PHC_1200ms_3.nii.gz
sub-01_ses-1_asl.nii.gz -> sub-01/ses-1/asl/sub-01_ses-1_asl.nii.gz
DEBUG: ep2d_se_pcasl_PHC_1200ms_3.nii.gz metadata: {'ASLContext': 'Label_Control', 'LabelingType': 'PCASL'}
DEBUG:
ep2d_se_pcasl_PHC_1200ms_MoCo_4.nii.gz
sub-01_ses-1_asl_moco.nii.gz -> sub-01/ses-1/asl/sub-01_ses-1_asl_moco.nii.gz
DEBUG:
bb11_faceA0_178_8.nii.gz
sub-01_ses-1_task-face_run-01_bold.nii.gz -> sub-01/ses-1/func/sub-01_ses-1_task-face_run-01_bold.nii.gz
DEBUG: bb11_faceA0_178_8.nii.gz metadata: {'TaskName': 'face'}
DEBUG:
bb11_faceB0_178_10.nii.gz
sub-01_ses-1_task-face_run-02_bold.nii.gz -> sub-01/ses-1/func/sub-01_ses-1_task-face_run-02_bold.nii.gz
DEBUG: bb11_faceB0_178_10.nii.gz metadata: {'TaskName': 'face'}
DEBUG:
bb11_cardB0_178_9.nii.gz
sub-01_ses-1_task-card_run-02_bold.nii.gz -> sub-01/ses-1/func/sub-01_ses-1_task-card_run-02_bold.nii.gz
DEBUG: bb11_cardB0_178_9.nii.gz metadata: {'TaskName': 'card'}
DEBUG:
ep2d_bold_MGH_13.nii.gz
sub-01_ses-1_task-rest_bold.nii.gz -> sub-01/ses-1/func/sub-01_ses-1_task-rest_bold.nii.gz
DEBUG: ep2d_bold_MGH_13.nii.gz metadata: {'TaskName': 'rest'}
(flywheel) T1nashes-Air:fw-heudiconv-demo mac3 [
```

# Walkthrough: Curating

6. Finally, when satisfied, remove the *dry\_run* flag, and apply your changes!

Don't be concerned if you mess up; **fw-heudiconv** only changes BIDS metadata and does not manipulate the underlying containers. Simply amend your heuristic and re-curate.

# Walkthrough: Curating

```
fw-heudiconv-demo — ttapera@chead:~ — -bash — 172x66
(flywheel) Tinashe-Air:fw-heudiconv-demo mac$ time fw-heudiconv-curate --project gear_testing --heuristic demo_heuristic.py --subject exampleSubject
/Users/mac/anaconda2/envs/flywheel/lib/python3.7/site-packages/fw_heudiconv/query.py:4: UserWarning: The DICOM readers are highly experimental, unstable, and only work for
Siemens time-series at the moment
Please use with caution. We would be grateful for your help in improving them
  from nibabel.nicom.dicomwrappers import wrapper_from_data
INFO: Querying Flywheel server...
INFO: Loading heuristic file...
INFO: Applying heuristic to query results...
Series not recognized: epi_singlerep epi_singlerep_5.nii.gz
Series not recognized: localizer_32channel localizer_32channel_1_i00003.nii.gz
Series not recognized: localizer_32channel localizer_32channel_1_i00002.nii.gz
Series not recognized: localizer_32channel localizer_32channel_1_i00001.nii.gz
Series not recognized: ep2d_effort1_236 ep2d_effort1_236_9.nii.gz
Series not recognized: ep2d_effort2_236 ep2d_effort2_236_10.nii.gz
Series not recognized: ep2d_effort3_1416 ep2d_effort3_1416_11.nii.gz
Series not recognized: ep2d_single ep2d_single_6.nii.gz
Series not recognized: localizer localizer_1_i00003.nii.gz
Series not recognized: localizer localizer_1_i00002.nii.gz
Series not recognized: localizer localizer_1_i00001.nii.gz
INFO: Processing IntendedFor fields based on heuristic file
INFO: Processing Metadata fields based on heuristic file
INFO: Applying changes to files...

real    0m16.875s
user    0m3.023s
sys     0m0.799s
(flywheel) Tinashe-Air:fw-heudiconv-demo mac$
```

# Walkthrough: Exporting

Optionally, use the export tool to view what your BIDS data looks like in a tree structure

```
$ fw-heudiconv-export --project gear_testing --path <path to where you want to download>  
--subject <optional list of subjects> --session <optional list of sessions> --dry_run <to  
print only a directory tree and not download>
```


# Walkthrough: Exporting

```
(flywheel) Tinashes-Air:fw-heudiconv-demo mac$ time fw-heudiconv-export --project gear_testing --path .
--subject exampleSubject --dry_run
/Users/mac/anaconda2/envs/flywheel/lib/python3.7/site-packages/fw_heudiconv/query.py:4: UserWarning: The
DICOM readers are highly experimental, unstable, and only work for Siemens time-series at the moment
Please use with caution. We would be grateful for your help in improving them
  from nibabel.nicom.dicomwrappers import wrapper_from_data
INFO: Gathering bids data:
INFO: Processing project files...
INFO: Processing session files...
INFO: Processing acquisition files...
INFO: Preparing output directory tree...
WARNING: No events.tsv found in func folder; creating empty TSVs
INFO: Done!
 bids_dataset/
  .bidsignore
  dataset_description.json
  sub-01/
    ses-2/
      anat/
        sub-01_ses-2_T1w.json
        sub-01_ses-2_T1w.nii.gz
    ses-1/
      anat/
        sub-01_ses-1_T1w.json
        sub-01_ses-1_T1w.nii.gz
      func/
        sub-01_ses-1_task-face_run-02_events.tsv
        sub-01_ses-1_task-card_run-02_events.tsv
        sub-01_ses-1_task-face_run-01_bold.json
        sub-01_ses-1_task-rest_bold.json
        sub-01_ses-1_task-card_run-02_bold.json
        sub-01_ses-1_task-face_run-02_bold.json
        sub-01_ses-1_task-face_run-02_bold.nii.gz
        sub-01_ses-1_task-card_run-02_bold.nii.gz
        sub-01_ses-1_task-rest_bold.nii.gz
        sub-01_ses-1_task-face_run-01_events.tsv
        sub-01_ses-1_task-face_run-01_bold.nii.gz

real    0m8.494s
user    0m1.972s
sys     0m0.760s
(flywheel) Tinashes-Air:fw-heudiconv-demo mac$
```

# Overview

1. Flywheel + BIDS = Awesome reproducible neuroscience!
2. Use **fw-heudiconv** to get your data BIDS ready
3. Community engagement is needed; submit **issues**, **bugs**, and **questions** to **Github**; **pull requests** welcome!
4. Show off your heuristics tips and tricks in our Repo too!



[github.com/PennBBL/fw-heudiconv/](https://github.com/PennBBL/fw-heudiconv/)



# Interactive workshop

Have a project on Flywheel you want  
to curate?

Stick around and we'll help you get  
started!



# Thanks to...

Ted Satterthwaite

Matt Cieslak

Azeez Adebimpe

Anna Xu

Diego Davila

Brain Behavior Lab

Desmond Oathes and Lab members

Matt Flounders

Joe Deluisi

UPenn Flywheel users

Flywheel support team